



Sandy L. Burkart PT, OCS, PhD is an orthopedic certified specialist recognized by the American Physical Therapy Association's Section of Orthopedic Physical Therapy for his teaching excellence. He presents evidenced based clinical workshops for physical and occupational therapists that integrate movement science with orthopedic examination algorithms, diagnostic specific

rehabilitation protocols, manual therapy practice, joint stabilization exercises and task specific training exercises that are functionally incorporated into the patient's activities of daily living. For workshop information go to: [www.functionalrehab.com](http://www.functionalrehab.com)



optp exclusive



*"The UE RANGER™ is one of the most innovative rehabilitation tools a therapist can use to improve functional range-of-motion, strength and coordination of both the orthopedic or neurologically impaired shoulder and upper extremity. The UE RANGER™ is ideal for post surgical shoulder rehabilitation or any acute shoulder dysfunction as it can be used for passive, active assistive, active and resistive range-of-motion in functional planes at varying degrees of elevation, rotation and scapular orientation."*

—Sandy L. Burkart, PhD, PT  
Functional Rehabilitation Seminars

**UE Ranger™**

The UE (Upper Extremity) Ranger™ is an assistive motion device designed to compliment the natural physiological process and movement health relations of the full upper extremity. The UE Ranger™ provides a progressive therapeutic balance between the necessary support and challenge requirements of patients striving to efficiently return to their full health and functional capacities. Features a molded plastic bilateral hand support and securing strap, telescopic tubing, articulating joints and detachable base. Included is a 30 page illustrated Clinical Manual and Home Exercise Program.

Ideal for:

- Passive assisted ROM performed by patient (never before possible)
- Progressive active assisted ROM
- Open and closed chain activity in all planes
- Functional neuromuscular retraining

UE Ranger™ (#435)

Video clips available at [www.optp.com](http://www.optp.com)



**UE Ranger™ Wall Mount**

Designed to facilitate and challenge a patient's capacity to produce functional movements using closed kinetic chain techniques. The stylish UE Ranger™ Wall Mount tracking frame easily mounts within your clinic. The base of the UE Ranger™ conveniently adapts and secures to the wall mount for a smooth, effortless transition during patients' treatment sessions. Non-returnable.

UE Ranger™ Wall Mount (#436)

**Upper Extremity Movement Health System**

The Upper Extremity Movement Health System comes complete with the UE Ranger™, Wall Mount, Clinical Manual and Home Exercise Program. Non-returnable.

Upper Extremity Movement Health System (#437PKG)



**UE Ranger™ Door Mount**

The door mount offers the benefits of the clinical wall mount, with:

- Versatility of installation
- Application with your patients Home Exercise Program

UE Ranger™ Door Mount (#438)

